

Human RBBP4 Recombinant Protein

Catalog Number: PROTQ09028

Overview

Product Name	Human RBBP4 Recombinant Protein
Description	Human RBBP4 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-425aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol, 1mM EDTA

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q09028, 1-425aa; MGSSHHHHHH SGLVPRGSH MGSMDADKEA AFDDAVEERV INEEYKIWKK NTPFLYDLVM THALEWPSLT AQWLDPVTRP EGKDFSIHRL VLGHTSDEQ NHLVIASVQL PNDDAQFDAS HYDSEKGEFG GFGSVSGKIE IEIKINHEGE VNRARYMPQN PCIIATKTPS SDVLVFDYTK HPSKPDPSGE CNPDLRLRGH QKEGYGLSWN PNLSGHLLSA SDDHTICLWD ISAVPKEGKV VDAKTIFTGH TAVVEDVSWH LLHESLFGSV ADDQKLMIWD TRSNNTSKPS HSDAHTAEV NCLSFNPYSE FILATGSADK TVALWDLRNL KKLHFSFESH KDEIFQVQWS PHNETILASS GTDRRLNVWD LSKIGEEQSP EDAEDGPPEL LFIHGHTAK ISDFSWNPNE PWWICSVSED NIMQVWQMAE NIYNDEDPEG SVDPEGQGS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human RBBP4 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.