

Human recombinant Hyaluronidase 1/HYAL1 protein

Catalog Number: PROTQ12794-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant Hyaluronidase 1/HYAL1 protein
Description	Human recombinant Hyaluronidase 1/HYAL1 protein expressed in HEK293 with His-Tag. Sequence domain: 22-435aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol.

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q12794, 22-435aa; FRGPLLPNRP FTTVWNANTQ WCLERHGVDV DVSVFDVVAN PGQTFRGPD M TIFYSSQLGT YPYTPTGEP VFGGLPQNAS LIAHLARTFQ DILAAIPAPD FSLAVIDWE AWRPRWAFNW DTKDIYRQRS RALVQAQHPD WPAPQVEAVA QDQFQGAARA WMAGTLQLGR ALRPRGLWGF YGFPCDNYD FLSPNYTGQC PSGIRAQNDQ LGWLWGQSRA LPSIYMPAV LEGTGKSQMY VQHRVAEAFR VAVAAGDPNL PVLPIYQIFY DTTNHFLPLD ELEHSLGES AQAAGVVLW VSWENTRTKE SCQAIKEYMD TTLGPFILNV TSGALLCSQA

LCSGHGRCVR RTSHPKALL LNPASFSIQL TPGGGPLSLR GALSLEDQAQ MAVEFKCRCY PGWQAPWCER KSMW

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant Hyaluronidase 1/HYAL1 protein](#)

For Research Use Only. Not for use in diagnostic procedures.