

## Human ELAVL2 Recombinant Protein

Catalog Number: PROTQ12926

### Overview

Product Name	Human ELAVL2 Recombinant Protein
Description	Human ELAVL2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-346aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 30% glycerol, 1mM DTT

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, Q12926, 1-346aa; MGSSHHHHHH SGLVPRGSH MGSMEQLSN GPTCNNTANG PTTINNNCSS PVDSGNTEDS KTNLIVNYLP QNMTQEELKS LFGSIGEIES CKLVRDKITG QSLGYGFVNY IDPKDAEKAI NTLNGLRLQT KTIKVSYPARP SSASIRDANL YVSGLPKMTMT QKELEQLFSQ YGRIITSRIL VDQVTGISRG VGFIRFDKRI EAEEAIKGLN GQKPPGATEP ITVKFANNPS QKTNQAILSQ LYQSPNRRYP GPLAQQQRF RFSPTMIDGM TSLAGINIPG HPGTGWCIFV YNLAPDADES ILWQMFGPFG AVTNVKVIRD FNTNKCKGFG FVTMTNYDEA AMAIASLNGY RLGDRVLQVS FKTNKTHKA

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human ELAVL2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.