

MRPL28 (NM_006428) Human Recombinant Protein

Catalog Number: PROTQ13084

Overview

Product Name	MRPL28 (NM_006428) Human Recombinant Protein
Description	Recombinant protein of human mitochondrial ribosomal protein L28 (MRPL28), nuclear gene encoding mitochondrial protein
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MPLHKYPVWLWKRLQLREGICSR L P G H Y L R S L E E E R T P T P V H Y R P H G A K F K I N P K N G Q R R E R V D V P I P I Y F P P E S Q R G L W G G E G W I L G Q I Y A N N D K L S K
R L K K V W K P Q L F E R E F Y S E I L D K K F T V T V M R T L D L I D E A Y G L D F Y I L K T P K E D L C S K F G M D L K R G M L L R L A R Q D P Q L H P E D P E R R A I Y D K Y K E F A I P E E E
A E W V G L T L E E A I E K Q R L L E E K D P V P L F K I Y V A E L I Q Q L Q Q A L S E P A V V Q K R A S G Q

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



MRPL28 (NM_006428) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.