

PAFAH (PLA2G7) (NM_001168357) Human Recombinant Protein

Catalog Number: PROTQ13093

Overview

Product Name	PAFAH (PLA2G7) (NM_001168357) Human Recombinant Protein
Description	Purified recombinant protein of human phospholipase A2, group VII (platelet-activating factor acetylhydrolase, plasma) (PLA2G7 / Lp-PLA2), transcript variant 2, with C-terminal DDK/His tag, expressed in human cells, 20 µg
Size	20 µg
Tag	C-DDK/His
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	PBS, pH7.4, 10% glycerol

Concentration

>50 µg/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

A DNA sequence encoding the region (Phe22 - Asn441) of human PLA2G7

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



PAFAH (PLA2G7) (NM_001168357) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.