

Human SLP-76/LCP2 Recombinant Protein

Catalog Number: PROTQ13094

Overview

Product Name	Human SLP-76/LCP2 Recombinant Protein
Description	Human SLP-76/LCP2 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 1-533aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q13094, 1-533aa; MALRNVP FRSEVLGWDP DSLADYFKKL NYKDCEKAVK KYHIDGARFL NLTENDIQKF PKLRVPILSK LSQEINKNEE
RRSIFTRKPV VPRFPEETES HEEDNGGWSS FEEDDYESPN DDQDGEDDGD YESPNEEEEA PVEDDADYEP PPSNDEEALQ NSILPAKPF
NSNSMYIDRP PSGKTPQPP VPPQRPMAAL PPPPAGRHS PLPPPQTNHE EPSRSRNHKT AKLPAPSIDR STKPPLDRSL APFDREPFTL
GKKPPFSDKP SIPAGRSLGE HLPKIQKPL PPTTERHERS SPLPGKKPPV PKHGWGPDRR ENEDDDVHQR PLPQPALLPM SSNTFPSRST
KPSPMNPLPS SHMPGAFSES NSSFPQSASL PPYFSQGPSN RPIRAEGRN FPLPLPNKPR PPSPAEEENS LNEEWYVSYI TRPEAEAALR
KINQDGTFLV RDSKKTITN PYVLMVLYKD KVYNIQIRYQ KESQVYLLGT GLRGKEDFLS VSDIIDYFRK MPLLLIDGKN RGSRYQCTLT HAAGYP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human SLP-76/LCP2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.