

Human TDP-43/TARDBP Recombinant Protein

Catalog Number: PROTQ13148

Overview

Product Name	Human TDP-43/TARDBP Recombinant Protein
Description	Human TDP-43/TARDBP Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-260aa. Application(s): SDS-PAGE.
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol 0.1M NaCl

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q13148, 1-260aa; MSEY IRVTEDEDE PIEIPSEDDG TVLLSTVTAQ FPGACGLRYR NPVSQCMRGV RLVEGILHAP DAGWGNLVVYV VNYPKDNKRK MDETDAASAV KVKRAVQKTS DLIVLGLPWK TTEQDLKEYF STFGEVLMVQ VKKDLKTGHS KGFGFVRFTE YETQVKVMSQ RHMIDGRWCD CKLPNSKQSQ DEPLRSRKVF VGRCTEDMTE DELREFFSQY GDVMDVFIPK PFRAFAVTF ADDQIAQSLC GEDLIKGIS VHSNA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human TDP-43/TARDBP Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.