

Human Cyclophilin G/PPIG Recombinant Protein

Catalog Number: PROTQ13427

Overview

Product Name	Human Cyclophilin G/PPIG Recombinant Protein
Description	Human Cyclophilin G/PPIG Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-175aa. Application(s): SDS-PAGE, Enzyme Activity.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q13427, 1-175aa; MGSSHHHHHH SSGLVPRGSH MGIKVRPRC FFDIAINNQP AGRVVFELFS DVCPKTCENF RCLCTGEKGT
GKSTQKPLHY KSCLFHRVVK DFMVQGGDFS EGNGRGGESI YGGFFEDES AVKHNKEFL SMANRGKDTN GSQFFITTKP TPHLDGHHVV
FGQVISGQEV VREIENQKTD AASKPFAEVR ILSCG

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Cyclophilin G/PPIG Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.