

Human NFATC2 Recombinant Protein

Catalog Number: PROTQ13469

Overview

Product Name	Human NFATC2 Recombinant Protein
Description	Human NFATC2 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 396-678aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 40% glycerol, 1mM DTT

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q13469, 396-678aa; PLEWPLSSQS GSYELRIEVQ PKPHHRAHYE TEGSRGAVKA PTGGHPVVQL HGYMENKPLG LQIFIGTADE RILKPHAFYQ VHRITGKTVT TTSYEKIVGN TKVLEIPLP KNNMRATIDC AGILKLRNAD IELRKGETDI GRKNTRVRLV FRVHIPESG RIVSLQTASN PIECSQRSAH ELP MVERQDT DSCLVYGGQQ MILTGQNFTS ESKVVFTEKT TDGQQIWEME ATVDKDKSQP NMLFVEIPEY RNKHIRTPVK VNFYVINGKR KRSQPQHFTY HPV

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human NFATC2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.