

Human IL-18 R alpha /IL-1 R5 Recombinant Protein

Catalog Number: PROTQ13478

Overview

Product Name	Human IL-18 R alpha /IL-1 R5 Recombinant Protein
Description	Human IL-18 R alpha /IL-1 R5 Recombinant Protein expressed in Baculovirus with hIgG-His-tag. Sequence domain: 19-329aa. Application(s): SDS-PAGE. Endotoxin: < 0.1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q13478, 19-329aa; AESCTSRPHI TVVEGEPFYL KHCSCLAHE IETTTKSWYK SSGSQEHVEL NPRSSRIAL HDCVLEFWPV ELNDTGSYFF QMKNYQKWK LNVIRRNKHS CFTERQVTSK IVEVKKFFQI TCENSYYQTL VNSTSLYKNC KKLLLENNKN PTIKKNAEFE DQGYYSVCVHF LHHNGKLFNI TKTFNITIVE DRSNIVPVLL GPKLNVHVAVE LGKNVRLNCS ALLNEEDVIY WMFGEENGSD PNIHEEKEMR IMTPEGKWHA SKVLRNIENIG ESNLNVLYNC TVASTGGTDT KSFILVRKAD MADIPGHVFT RVEPKSCDKT HTCPPCPAPE LLGGPSVFLF PPKPKDTLMI SRTPEVTCVV VDVSHEDPEV KFNWYVDGVE VHNAKTKPRE EQYNSTYRVV SVLTVLHQDW LNGKEYKCKV SNKALPAPIE KTISKAKGQP REPQVYTLPP SRDELTKNQV SLTCLVKGFY PSDIAVEWES NGQPENNYKT TPPVLDSDGS FFLYSKLTVD KSRWQQGNVF SCSVMHEALH NHYTKSLSL SPGKHHHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human IL-18 R alpha /IL-1 R5 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.