

Human Smad4 Recombinant Protein

Catalog Number: PROTQ13485

Overview

Product Name	Human Smad4 Recombinant Protein
Description	Human Smad4 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-552aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q13485, 1-552aa; MGSSHHHHHH SGLVPRGSH MDNMSITNTP TSNDACLSIV HSLMCHRQGG ESETFAKRAI ESLVKKLKEK
KDELDSLITA ITTNGAHPK CVTIQRTLDG RLQVAGRKGF PHVIYARLWR WPD LHKNELK HVKYCQYAFD LKCD SVCVNP YHYERVVSPG
IDLSGLTLQS NAPSSMMVKD EYVHDFEGQP SLSTEGHSIQ TIQHPPSNRA STETYSTPAL LAPSES NATS TANFPNIPVA STSQPASILG GSHSEGLLQI
ASGPQPGQQQ NGFTGQPATY HHNSTTTWTG SRTAPYTPNL PHHQNGHLQH HPPMPHPGH YWPVHNELAF QPPISNHPAP EYWCSIAYFE
MDVQVGETFK VPSSCPIVTV DGYVDPSGGD RFCLGQLSNV HRTEAIERAR LHIGKGVQLE CKGEDVWVWR CLSDHAVFVQ SYLDREAGR
APGDAVHKIY PSAYIKVFDL RQCHRQMQQQ AATAQAAAAA QAAAVAGNIP GPGSVGGIAP AISLSAAAGI GVDDLRLRCI LRMSFVKGWG
PDYPRQSIKE TPCWIEIHLH RALQLLDEVL HTMPIADPQP LD

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Smad4 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.