

Human Dynactin Subunit 2/DCTN2 Recombinant Protein

Catalog Number: PROTQ13561

Overview

Product Name	Human Dynactin Subunit 2/DCTN2 Recombinant Protein
Description	Human Dynactin Subunit 2/DCTN2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-401aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.15M NaCl, 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q13561, 1-401aa; MGSSHHHHHH SGLVPRGSH MGSMADPKYA DLPGIARNEP DVYETSDLPE DDQAEFDAEE LTSTSVEHII
VNPNAAYDKF KDKRVGKGL DFSDRIGKTK RTGYESGEYE MLGEGLVKVE TPQQKYQRLL HEVQELTTEV EKIKTTVKES ATEEKLTPVL
LAKQLAALKQ QLVASHLEKL LGPDAAINLT DPDGALAKRL LLQLEATKNS KGGSGGKTTG TPPDSSLVTY ELHSRPEQDK FSQAAKVAEL
EKRLTELETA VRCDQDAQNP LSAGLQGAQL METVELLQAK VSALDLAVLD QVEARLQSVL GKVNEIAKHK ASVEDADTQS KVHQLYETIQ
RWSPIASTLP ELVQRLVTIK QLHEQAMQFG QLLTHLDTTQ QMIANSLKDN TLLTQVQTT MRENLATVEG NFASIDERMK KLGK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Dynactin Subunit 2/DCTN2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.