

Human CuL1 Recombinant Protein

Catalog Number: PROTQ13616

Overview

Product Name	Human CuL1 Recombinant Protein
Description	Human CuL1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-410aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol, 100mM NaCl

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q13616, 1-410aa; MSSTRSQNPH GLKQIGLDQI WDDLRAIQQ VYTRQSMAS RYMELYTHVY NYCTSVHQSN QARGAGVPPS
KSKKGQTPGG AQFVGLLELYK RLKEFLKKNYL TNLLKDGEDL MDESVLKFYT QQWEDYRFSS KVLNGICAYL NRHWVRRECD EGRKGIYEIY
SLALVTWRDC LFRPLNKQVT NAVLKIEKE RNGETINTRL ISGVVQSYVE LGLNEDDAFA KGPTLVYKE SFESQFLADT ERFYTRESTE
FLQQNPVTEY MKKAEARLLE EQRRVQVYLH ESTQDELARK CEQVLEKHL EIFHTEFQNL LDADKNE DLG RMYNLVSRIQ DGLGELKLL
ETHIHNQGLA AIEKCGEAL NDPKMYVQTV LDVHKYNAL VMSAFNNDAG FVAALDKACG RFINNNAVTK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CuL1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.