

## RAB31 (NM\_006868) Human Recombinant Protein

Catalog Number: PROTQ13636

### Overview

|              |  |
|--------------|--|
| Product Name | RAB31 (NM_006868) Human Recombinant Protein                            |
| Description  | Recombinant protein of human RAB31, member RAS oncogene family (RAB31) |
| Size         | 20 µg  |
| Tag          | C-Myc/DDK  |
| Form         | Frozen Solution in PBS Buffer  |
| Source       | HEK293T  |
| Formulation  | 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol                   |

### Concentration

>50 ug/mL as determined by microplate BCA method

### Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

### Amino Acid Sequence

MMAIRELKVCLLGDTGVGKSSIVCRFVQDHFHDHNSPTIGASFMTKTVPCGNELHKFLIWDTAGQERFHS LAPMYRGSAAAIVYDITKQDSFYTLKK  
WVKELKEHGPENIVMAIAGNKCDLSDIREVPLKDAKEYAESIGAI VVETS AKNAINIEELFQGISRQIPPLDPHENGNGTIKVEKPTMQASRRCC

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



RAB31 (NM\_006868) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.