

## Human ALCAM/CD166 Recombinant Protein

Catalog Number: PROTQ13740

### Overview

Product Name	Human ALCAM/CD166 Recombinant Protein
Description	Human ALCAM/CD166 Recombinant Protein expressed in Baculovirus with hlgG-His-tag. Sequence domain: 28-527aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, Q13740, 28-527aa; ADPWYTVNSA YGDTIIIPCR LDVPQNL MFG KWKYEKPDGS PVFIAFRSST KKS VQYDDVP EYKDRNLNLS  
NYTLSISNAR ISDEKRFVCM LVTEDNVFEA PTIVKVKQP SKPEIVSKAL FLETEQLKKL GDCISEDSYP DGNITWYRNG KVLHPLEGAV VIIFKKEMDP  
VTQLYTMTST LEYKTTKADI QMPFTCSVTY YGPGSQKTIH SEQAVFDIYY PTEQVTIQVL PPKNAIKEGD NITLKCLGNG NPPPEEFLFY  
LPGQPEGIRS SNTYTLTDVR RNATGDYKCS LIDKKS MIAS TAITVHYLDL SLNPSGEVTR QIGDALPVSC TISASRNATV VWMKDNIRLR SSPSFSSLHY  
QDAGNYVCET ALQEVEGLKK RESLTLIVEG KPQIKMTKKT DPSGLSKTII CHVEGF PKPA IQWTITGSGS VINQTEESPY INGRYYSKII ISPEENVTLT  
CTAENQLERT VNSLNV SANE NREKVNDQAK LIVGIVGLL LAALEPKSCD KTHTCP PCA PELLGGPSVF LFPPKPKDTL MISRTPEVTC  
VVVDVSHEDP EVKFNWYVDG VEVHNAKTKP REEQYNSTYR VVSVLTVLHQ DWLNGKEYKC KVS NKALPAP IEKTISKAKG QPREPQVYTL  
PPSRDELTKN QVSLTCLVKG FYPSDIAVEW ESNQGPENNY KTPPVLDSD GSFFLYSKLT VDKSRWQQGN VFSCSVMHEA LHNHYTQKSL  
SLSPGKHHHH HH

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human ALCAM/CD166 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.