

## Human recombinant ENPP-2/Autotaxin protein

Catalog Number: PROTQ13822

### Overview

Product Name	Human recombinant ENPP-2/Autotaxin protein
Description	Human recombinant ENPP-2/Autotaxin protein expressed in HEK293 with His-Tag. Sequence domain: 49-863aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.25mg/ml (determined by Absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, Q13822, 49-863aa; MDSPWTN ISGCKGRFC ELQEAGPPDC RCDNLCKSYT SCCHDFDEL C LKTARGWECT KDRCGEVRNE ENACHCSEDC LARGDCCTNY QVVKGESHW VDDDC EEIKA AEC PAGFVRP PLIIFSVDGF RASYM KKGSK VMPNIEKLRS CGTHSPYMRP VYPTKTFPNL YTLATGLYPE SHGIVGNSMY DPVFDATFHL RGREKFNHRW WGGQPLWITA TKQGVKAGTF FWSVVIPHER RILTILQWLT LPDHERPSVY AFYSEQPDFS GHKYGPFGE MTNPLREIDK IVGQLMDGLK QKLRHRCVNV IFVGDHGMED VTCDRTEFLS NYLTNVDDIT LVPGLGRIR SKFSNNAKYD PKAIANLTC KKPQDHF KP Y LKQHLPKRLH YANNRRIEDI HLLVERRWHV ARKPLDVYKK PSGKCFQGD HGF DNKVN SM QTVFVG YGST FKYKTKVPPF ENIELYNVMC DLLGLKPAPN NGTHGSLNHL LRTNTFRPTM PEEVTRPNYP GIMYLQSDFD LGCTCDDKVE PKNKLDL NK RLHTKGSTEE RHLLYGRPAV LYRTRYDILY HTDFESGYSE IFLMPLWTSY TVSKQAEVSS VPDHLTSCVR PDVRVSPSFS QNCLAYKNDK QMSYGFLFPP YLSSSPEAKY DAFLVTNMVP MYP AFKRVWN YFQRVLVKKY ASERNGV NVI SGPIFDYDYD GLHDTE DKIK QYVEGSSIPV PTHYYSIITS CLDFTQPADK CDGPLSVSSF ILPHRPDNEE SCNSEDESK WVEELMKMHT ARVRDIEHLT SLDFFRKT SR SYPEILTLKT YLHTYESEI

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant ENPP-2/Autotaxin protein

For Research Use Only. Not for use in diagnostic procedures.