

Human GPNMB Recombinant Protein

Catalog Number: PROTQ14956

Overview

Product Name	Human GPNMB Recombinant Protein
Description	Human GPNMB Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 22-474aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by BCA assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q14956, 22-474aa; AKRFHDV LGNERPSAYM REHNQLNGWS SDENDWNEKL YPVWKRGD MR WKNSWKGG RV QAVLTS DSPA LVGSNITFAV NLIFPRCQKE DANGNIVYEK NCRNEAGLSA DPVYVNWTA W SEDSDGENGT GQSHHNVFPD GKPFPHHPGW RRWNFIYVFH TLGQYFQKLG RCSVRVSVNT ANVTLGPQLM EVT VYRRHGR AYVPIAQVKD VYVVT DQIPV FVTMFQKNDR NSSDETFLKD LPIMFDVLIH DPSHFLNYST INYKWSFGDN TGLFVSTNHT VNHTYVLNGT FSLNLT VKAA APGPCPPPPP PPRPSKPTPS LGPAGDNPLE LSRIPDENCQ INRYGHFQAT ITIVEGILEV NIIQMTDVL M PVPWPESLI DFVVT CQGSI PTEVCTIISD PTCEITQNTV CSPVDVDEMC LLTVRRTFNG SGT YCVNLTL GDDTSLALTS TLISVP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human GPNMB Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.