

## Human BRD3 Recombinant Protein

Catalog Number: PROTQ15059

### Overview

Product Name	Human BRD3 Recombinant Protein
Description	Human BRD3 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-416aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, Q15059, 1-416aa; MSTATTV APAGIPATPG PVNPPPPEVS NPSKPRKTN QLQYMQNVVV KTLWKHQFAW PFYQPVDAIK LNLDPDYHKII  
KNPMDMGTIK KRLENNYYWS ASECMQDFNT MFTNCYIYNK PTDDIVLMAQ ALEKIFLQKV AQMPQEEVEL LPPAPKGKGR KPAAGAQSAG  
TQQVAAVSSV SPATPFQSVP PTVSQTPIA ATPVPTITAN VTSVPVPPAA APPPPATPIV PVVPPTPPVV KKKGVKRRKAD TTTPTTSAIT ASRSESPPL  
SDPKQAKVVA RRESGGRPIK PPKKDLEDGE VPQHAGKKGK LSEHLRYCDS ILREMLSKKH AAYAWPFYKP VDAEALELHD YHDIKHPMD  
LSTVKKRMDG REYPDAQGFA ADVRLMFSNC YKYNPPDHEV VAMARKLQDV FEMRFAKMP

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human BRD3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.