

PD1 (PDCD1) (NM_005018) Human Recombinant Protein

Catalog Number: PROTQ15116

Overview

Product Name	PD1 (PDCD1) (NM_005018) Human Recombinant Protein
Description	Purified recombinant protein of Homo sapiens programmed cell death 1 (PD1/PDCD1), residues 25-167aa, with C-terminal Fc tag, expressed in HEK293 cells.
Size	20 µg
Tag	C-Fc
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	PBS, pH7.4, 10% glycerol

Concentration

>50 µg/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

A DNA sequence encoding the extracellular domain (Leu25 - Gln167) of human programmed cell death 1 (PD1)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



PD1 (PDCD1) (NM_005018) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.