

Human QPRT Recombinant Protein

Catalog Number: PROTQ15274

Overview

Product Name	Human QPRT Recombinant Protein
Description	Human QPRT Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-297aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q15274, 1-297aa; MGSSHHHHHH SGLVPRGSH MDAEGLALLL PPVTLAALVD SWLREDCPGL NYAALVSGAG PSQAALWAKS
PGVLGQPF DAIFTQLNCQ VSWFLPEGSK LVPVARVAEV RGAHCLLLG ERVALNTLAR CSGIASAAAA AVEAARGAGW TGHVAGTRKT
TPGFRLVEKY GLLVGAASH RYDLGGLVMV KDNHVVAAGG VEKAVRAARQ AADFALKVEV ECSSLQEAVQ AAEAGADLVL LDNFKPEELH
PTATVLKAQF PSVAVEASGG ITLDNLPQFC GPHIDVISMG MLTQAAPALD FSLKLFKEV APVPKIH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human QPRT Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.