

## Human QPRT Recombinant Protein

Catalog Number: PROTQ15274

### Overview

Product Name	Human QPRT Recombinant Protein
Description	Human QPRT Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-297aa. Application(s): SDS-PAGE.
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol

### Storage

Can be stored at +2°C to +8°C for 1 week. For long term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

MGSSHHHHHH SSGLVPRGSH MDAEGLALLL PPVTLAALVD SWLREDCPGL NYAALVSGAG PSQAALWAKS PGVLAGQPFF DAIFTQLNCQ  
VSWFLPEGSK LVPVARVAEV RGPACHLLLG ERVALNTLAR CSGIASAAAA AVEAARGAGW TGHVAGTRKT TPGFRLVEKY GLLVGGAASH  
RYDLGGLVMV KDNHVVAAGG VEKAVRAARQ AADFALKVEV ECSSLQEAQV AAEAGADLVL LDNFKPEELH PTATVLKAQF PSVAVEASGG  
ITLDNLPQFC GPHIDVISMG MLTQAAPALD FSLKLFKEV APVPKIH

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human QPRT Recombinant Protein