

Human recombinant ErbB4/Her4 protein

Catalog Number: PROTQ15303

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant ErbB4/Her4 protein
Description	Human recombinant ErbB4/Her4 protein expressed in HEK293 with hlgG-His-Tag. Sequence domain: 26-649aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q15303, 26-649aa; QSV CAGTENK LSSLDLEQQ YRALRKYYEN CEVVMGNLEI TSIEHNRDLS FLRSVREVTG YVLVALNQFR YLPLENLR II
RGTKLYEDRY ALAIFLN YRK DGNFGLQELG LKNLTEILNG GVVVDQNKFL CYADTIHWQD IVRNPWPSNL TLVSTNGSSG CGRCHK SCTG
RCWGP TENHC QTLTRTVCAE QCDGR CYGPY VSDCCHRECA GGCSGPKDTD CFACMNFNDS GACVTQCPQT FVYNPTTFQL EHNFN AKYTY
GAFCVKKCPH NFVVDSSSCV RACPSSKMEV EENG IKMCKP CTDICPKACD GIGTGLMSA QTV DSSNIDK FINCTKINGN LIFLVTGIHG

DPYNAIEAID PEKLVNVRTV REITGFLNIQ SWPPNMTDFS VFSNLVTIGG RVLYSGLSLL ILKQQGITSL QFQSLKEISA GNIYITDNSN LCYYHTINWT
TLFSTINQRI VIRDNRKAEN CTAEGMVCNH LCSSDGCWGP GPDQCLSCRR FSRGRICIES CNLYDGEFRE FENGSI CVEC DPQCEK MEDG
LLTCHGPGPD NCTKCSHFKD GPNCVEKCPD GLQGANSFIF KYADPDRECH PCHPNCTQGC NGPTSHDCIY YPWTGHSTLP QHAR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant ErbB4/Her4 protein](#)

For Research Use Only. Not for use in diagnostic procedures.