

Human RheB Recombinant Protein

Catalog Number: PROTQ15382

Overview

Product Name	Human RheB Recombinant Protein
Description	Human RheB Recombinant Protein expressed in E. coli with T7-tag. Sequence domain: 1-181aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	T7-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q15382, 1-181aa; MASMTGGQQM GRGSASMPQS KSRKIALGY RSVGKSSLTI QFVEGQFVDS YDPTIENTFT KLITVNGQEY HLQLVDTAGQ DEYSIFPQTY SIDINGYILV YSVTSIKSFE VIKVIHGKLL DMVGKVQIPI MLVGNKKDLH MERVISYEEG KALAESWNAA FLESSAKENQ TAVDVFRRIL LEAEKMDGAA SQGKSSC

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human RheB Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.