

Human beta IG-H3 Recombinant Protein

Catalog Number: PROTQ15582

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human beta IG-H3 Recombinant Protein
Description	Human beta IG-H3 Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 502-636aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0)

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q15582, 502-636aa; MGTVM DVLKG DNRFSMLVAA IQSAGLTETL NREGVYTVFA PTNEAFRALP PRERSRL LGD AKELANILKY
HIGDEILVSG GIGALVRLKS LQGD KLEVSL KNNVSVNKE PVAEPDIMAT NGVVHVITNV LQPPA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human beta IG-H3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.