

## Human MAPRE1 Recombinant Protein

Catalog Number: PROTQ15691

### Overview

Product Name	Human MAPRE1 Recombinant Protein
Description	Human MAPRE1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-268aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, Q15691, 1-268aa; MAVNVYSTSV TSDNLSRHDM LAWINESLQL NLTKIEQLCS GAAYCQFMDM LFGSIALKK VKFQAKLEHE YIQNFKILQA GFKRMGVDKI IPVDKLVKGG FQDNFEFVQW FKKFFDANYD GKDYDPVAAR QGQETAVAPS LVAPALNPKK KPLTSSSAAP QRPISQRTA AAPKAGPGVV RKNPGVGNPD DEAAELMQQV NVLKLTVEDL EKERDFYFGK LRNIELICQE NEGENDPVLQ RIVDILYATD EGFVIPDEGG PQEEQEEY [?/description\\_after\\_attributes?](#)

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human MAPRE1 Recombinant Protein