

Human NSDHL Recombinant Protein

Catalog Number: PROTQ15738

Overview

Product Name	Human NSDHL Recombinant Protein
Description	Human NSDHL Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-297aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol 0.1M NaCl

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 80% by SDS-PAGE

Amino Acid Sequence

Human, Q15738, 1-297aa; MGSSHHHHHH SSGLVPRGSH MGSMEPAVSE PMRDQVARTH LTEDTPKVNA DIEKVNQNQA KRCTVIGGSG
FLGQHMVEQL LARGYAVNVF DIQQGFDNPQ VRFFLDLCS RQDLYPALKG VNTVFHCASP PPSSNNKELF YRVNYIGTKN VIETCKEAGV
QKLILTSSAS VIFEGVDIKN GTEDLPYAMK PIDYYTETKI LQERAVLGAN DPEKNFLT TA IRPHGIFGPR DPQLVPILIE AARNGKMKFV IGNGKNLVDF
TFVENVVHGH ILAAEQLSRD STLGGKAFHI TNDPIPFWT FLSRILTGLN YEAPKYHIPY

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human NSDHL Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.