

Human NEDD8 Recombinant Protein

Catalog Number: PROTQ15843

Overview

Product Name	Human NEDD8 Recombinant Protein
Description	Human NEDD8 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-76aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.5mM DTT, 50mM NaCl, 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q15843, 1-76aa; MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSHMLI KVKTLTGKEI EIDIEPTDKV ERIKERVEEK
EGIPPQQRL IYSGKQMND E KTAADYKILG GSVLHLVLAL RGG ?/description_after_attributes?

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human NEDD8 Recombinant Protein