

Human IGFBP-rp1/IGFBP-7 Recombinant Protein

Catalog Number: PROTQ16270

Overview

Product Name	Human IGFBP-rp1/IGFBP-7 Recombinant Protein
Description	Human IGFBP-rp1/IGFBP-7 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 27-282aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% glycerol, 2mM DTT, 1mM EDTA

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q16270, 27-282aa; MGSSHHHHHH SGLVPRGSH MGSSSDTCG PCEPASCPL PPLGCLLGET RDACGCCPMC ARGEGEPCGG
GGAGRGYCAP GMECVKSRKR RKGKAGAAAG GPGVSGVCVC KSRYPVCGSD GTTYPGCGQL RAASQRAESR GEKAITQVSK GTCEQGPSIV
TPPKDIWNVT GAQVYLSCEV IGIPTPLIWNKVKRGRHYGV QRTELLPGDR DNLAIQTRGG PEKHEVTGWV LVSPLSKEDA GEYECHASNS
QGQASASAKI TVVDALHEIP VKKGEGAEI

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human IGFBP-rp1/IGFBP-7 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.