

Human SGCB Recombinant Protein

Catalog Number: PROTQ16585

Overview

Product Name	Human SGCB Recombinant Protein
Description	Human SGCB Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 87-318aa. Application(s): SDS-PAGE, Denatured.
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Storage

Can be stored at +2°C to +8°C for 1 week. For long term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

WAVIRIG PNGCDSMEFH ESGLLRFKQV SDMGVIHPLY KSTVGGRRNE NLVITGNNQP IVFQQGTTKL SVENNKTSIT SDIGMQFFDP RTQNILFSTD
YETHEFHLPS GVKSLNVQKA STERITSNAT SDLNIKVDGR AIVRGNEGVF IMGKTIEFHM GGNMELKAEN SIILNGSVMV STTRLPSSSS
GDQLGSGDWV RYKLCMCADG TLFKVQVTSQ NMGCQISDNP CGNTH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human SGCB Recombinant Protein