

## Human Frataxin Recombinant Protein

Catalog Number: PROTQ16595

### Overview

|              |  |
|--------------|--|
| Product Name | Human Frataxin Recombinant Protein   |
| Description  | Human Frataxin Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 42-210aa. Application(s): SDS-PAGE. |
| Size         | Starting from 20ug   |
| Tag          | His-Tag  |
| Form         | Liquid   |
| Source       | E. coli  |
| Formulation  | 20mM Tris-HCl buffer (pH 8.5) containing 20% glycerol, 0.1M NaCl 1mM DTT   |

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Human, Q16595, 42-210aa; LRTDIDATC TPRRASSNQR GLNQIWNVKK QSVYLMNLRK SGTLGHPGSL DETTYERLAE ETLDSLAEFF EDLADKPYTF EDYDVSGSG VLTVKLGGDL GTYVINKQTP NKQIWLSSPS SGPKRYDWTG KNWVYSHDGV SLHELLAAEL TKALKTKLDL SSLAYSGKDA

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human Frataxin Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.