

Human recombinant TrkB protein

Catalog Number: PROTQ16620-2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant TrkB protein
Description	Human recombinant TrkB protein expressed in HEK293 with hIgG-Tag. Sequence domain: 32-430aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

1mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q16620, 32-430aa; CPTSCCKSAS RIWCDSPSPG IVAFPRLEPN SVDPENITEI FIANQKRLEI INEDDVEAYV GLRNLTI VDS GLKFVAHKAF LKNSNLQHIN FTRNKLTSL RKHFRHLDLS ELILVGNPFT CSCDIMWIKT LQEA KSSPD T QDLYCLNESS KNIPLANLQI PNCGLPSANL AAPNLTV EEG KSI TLSCSVA GDPV PNM YWD VGNLVSKHMN ETSHTQGSLR ITNISSDDSG KQISCVAENL VGEDQDSVNL TVHFAPTITF LESPTSDHHW CIPFTVKG NP KPALQWFYNG AILNESKYIC TKIHVTNHT E YHGCLQLDNP THMNNGDYTL IAKNEYGKDE KQISAHFMGW PGIDDGANPN

YPDVIYEDYG TAANDIGDTT NRSNEIPSTD VTDKTGREH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant TrkB protein](#)

For Research Use Only. Not for use in diagnostic procedures.