

Human Syntaxin 1A Recombinant Protein

Catalog Number: PROTQ16623

Overview

Product Name	Human Syntaxin 1A Recombinant Protein
Description	Human Syntaxin 1A Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-265aa. Application(s): SDS-PAGE.
Size	Starting from 4ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q16623, 1-265aa; MKDRTQELRT AKDSDDDDDV AVTVDRDRFM DEFFEQVEEI RGFIDKIAEN VEEVKRKHSA ILASPNPDEK
TKEELEELMS DIKKTANKVR SKLKSIEQSI EQEGLNRSS ADLRIRKTQH STLSRKFVEV MSEYNATQSD YRERCKGRIQ RQLEITGRRT TSEELEDMLE
SGNPAIFASG IIMDSSISKQ ALSEIETRHS EIIKLENSIR ELHDMFMDMA MLVESQGEMI DRIEYNVEHA VDYVERAVSD TTKAVKYQSK ARKK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Syntaxin 1A Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.