

Human Melan-A/MART-1 Recombinant Protein

Catalog Number: PROTQ16655

Overview

| | |
|--------------|--|
| Product Name | Human Melan-A/MART-1 Recombinant Protein |
| Description | Human Melan-A/MART-1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 48-118aa. Application(s): SDS-PAGE. |
| Size | Starting from 20ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT |

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q16655, 48-118aa; MGSSHHHHHH SSGLVPRGSH MGSCRRRNGY RALMDKSLHV GTQCALTRRC PQEGFDHRDS KVSLQEKNC E
PVVPNAPPAY EKLSAEQSP PYS P

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Melan-A/MART-1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.