

DEC1 (NM_001359) Human Recombinant Protein

Catalog Number: PROTQ16698

Overview

Product Name	DEC1 (NM_001359) Human Recombinant Protein
Description	Recombinant protein of human 2,4-dienoyl CoA reductase 1, mitochondrial (DEC1), nuclear gene encoding mitochondrial protein
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MKLPARVFFTLGSR L P C G L A P R R F F S Y G T K I L Y Q N T E A L Q S K F F S P L Q K A M L P P N S F Q G K V A F I T G G G T G L G K G M T T L L S S L G A Q C V I A S R K M D V L K A T A
E Q I S S Q T G N K V H A I Q C D V R D P D M V Q N T V S E L I K V A G H P N I V I N N A A G N F I S P T E R L S P N A W K T I T D I V L N G T A F V T L E I G K Q L I K A Q K G A A F L S I T T I Y A E T
G S G F V V P S A S A K A G V E A M S K S L A A E W G K Y G M R F N V I Q P G P I K T K G A F S R L D P T G T F E K E M I G R I P C G R L G T V E E L A N L A A F L C S D Y A S W I N G A V I K F D G
G E E V L I S G E F N D L R K V T K E Q W D T I E E L I R K T K G S

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



DEC1 (NM_001359) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.