

## CLPP (NM\_006012) Human Recombinant Protein

Catalog Number: PROTQ16740

### Overview

|              |  |
|--------------|--|
| Product Name | CLPP (NM_006012) Human Recombinant Protein   |
| Description  | Recombinant protein of human ClpP caseinolytic peptidase, ATP-dependent, proteolytic subunit homolog (E. coli) (CLPP), nuclear gene encoding mitochondrial protein |
| Source       | HEK293T  |
| Formulation  | 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol   |

### Storage

Store at -80°C.

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



CLPP (NM\_006012) Human Recombinant Protein