

Carbonic Anhydrase IX (CA9) (NM 001216) Human Recombinant Protein

Catalog Number: PROTQ16790

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Carbonic Anhydrase IX (CA9) (NM_001216) Human Recombinant Protein
Description	Purified recombinant protein of Human carbonic anhydrase IX (CA9), with C-terminal His tag, secretory expressed in HEK293 cells, 50ug
Size	50 µg
Tag	C-His
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	PBS, pH7.4, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

A DNA sequence encoding the region Gln38-Asp414 of CAIX

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Carbonic Anhydrase IX \(CA9\) \(NM_001216\) Human Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.