

Human MTHFSD Recombinant Protein

Catalog Number: PROTQ2M296

Overview

Product Name	Human MTHFSD Recombinant Protein
Description	Human MTHFSD Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-383aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 30% glycerol, 1mM DTT

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q2M296, 1-383aa; MGSSHHHHHH SSGLVPRGSH MGSMEPRAVG VSKQDIREQI WGYMESQNLA DFPRPVHHRI PNFKGASHAA
EQLPRLQAFK TARTIKVNPDPAPQKSARFFV LESKKTLLVP TPLRLTGLFN KITPPPATK DILRKCATSQ GVRNYSVPIG LDSRVLVDLV
VVGSAVSEK GWRIGKGEY ADLEYAMMVS MGAVSKETPV VTIVHDCQVV DIPEELVEEH DITVDYILTP TRVIATGCKR PKPMGITWFK
ISLEMMEKIP ILRSLRAREQ QAGKDVTLQG EHQHLPPEGC QQTVP LSVGR RPPDTPGPET NSMEAAPGSP PGEGAPLAAD VYVGNLPGDA
RVSDLKRALR ELGSVPLRLT WQGPRRRAFL HYPDSAAAQQ AVSCLQGLRL GTDTRLVALA RQQRDK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human MTHFSD Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.