

Mouse recombinant CD39L2/ENTPD6 Protein

Catalog Number: PROTQ3U0P5

Overview

Product Name	Mouse recombinant CD39L2/ENTPD6 Protein
Description	Mouse recombinant CD39L2/ENTPD6 Protein protein expressed in HEK293 with His-Tag. Sequence domain: 33-455aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, Q3U0P5, 33-455aa; KWHRAS AAQAFFTIAG AASGARWTQQ AFSSPGSAAR GHEVFGIMF DAGSTGTRIH VFQFARPPGE TPTLTHETFK ALKPGLSAYA DDVEKSAQGI QELLNVAKQH IPYDFWKATP LVLKATAGLR LLPGEKAQKL LQKVKEVFK SPFLVGDDCV SIMNGTDEGV SAWITVNFLT GSLKTPGSSS VGMLDLGGGS TQITFLPRVE GTLQASPPGH LTALQMFNRT YKLYSYSYLG LGLMSARLAI LGGVEGKPAE NDKELVSPCL SPRFRGEWEH AEVTYRISGQ KAVGLYELCA SRVSEVLRNK VHRTEEAQHV DFYAFSYYYD LAASFGLIDA EKGGSLLVVD FEIAAKYVCR TLETQPPSP FACMDLTYIS LLLHEFGFPG DKVLKLRKI DNVETSWALG AIFHYIDSLK RQKVPAL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse recombinant CD39L2/ENTPD6 Protein

For Research Use Only. Not for use in diagnostic procedures.