

TTC18 (CFAP70) (NM_145170) Human Recombinant Protein

Catalog Number: PROTQ5T0N1

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	TTC18 (CFAP70) (NM_145170) Human Recombinant Protein
Description	Recombinant protein of human tetratricopeptide repeat domain 18 (TTC18)
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MEQVPSAGRLVQITVTEGYDLKGFKGDTPVTFIRAEFNQVVLGDSAKITVSPEGSAKYNFTSSFEFNPEGGITSDDLAKHPVFLTVTEVLPKEKKQKEEKT
LILGQAVVDLLPLLEGQSSFQTTVPLHPVQGSPLTPRSSAKQCSLEVKVLVAEPLTTAQISGGNLLKVTLAAAYSVPESFIPTGPGQNYMVGVLQVPSLGE
KDYPILFKNGTLKLGGEREPVPRPKKWPIANILAPGANNIPDAFIVGGPYEEEEGELNHPEDSEFRNQAECIKKRIIWDLESRCYLDPSAVVSFQKRIADCR
LWPVEITRVPLVTIPKGKAGKTEKTDDEAQLSFHGVAYVNMVPLLYPGVKRIRGAFHVYPYLDSSVHEKTKCLLSLFRDIGHHLIHNNKIGGINSLLSKQA
VSKNLKEDKPVKEKDIDGRPRPGDVQAPSIKSQSSDTPLEGEPLSHNPEGQQYVEAGTYIVLEIQLDKALVPKRMPEELARRVKEMIPPRPLTRRTGGA
QKAVSDYHIQIKNISRAILDEYYRMFGKQVAKLESMDSETLEEQKCQLSYELNCSGKYFAFKEQLKHAVVKIVRD KYLKTTSFESQEELQTFISELYVFLV

DQMHVALNQTMPDDVQGTVATIYTSSEQLQLFAFEAEVNENFEMAAAYYKERLVREPQNLDHWLDYGAFCLLTEDNIKAQECFQKALSLNQSHIHSLL
LCGVLAVLLENYEQAEIFFEDATCLEPTNVVAWTLLGLYYEIQNNDIRMEMAFHEASKQLQARMLQAQVTKQKSTGVEDTEERGKRESSLGPWGITNGS
ATAIKVEAPAGPGAALSILDKFLEESSKLQSDSQEPILTTQTWDPSISQKPSNTFIKEIPTKKEASKCQDSSALLHPGLHYGVSQTTTIFMETIHFLMKVKAV
QYVHRVLAHELLCPQGGPSCEYYLVLAQTHILKKNFAKAEYLLQQAQMDYLNPNVWGLKGHLYFLSGNHSEAKACYERTISFVVDASEMHFIFLRLGL
IYLEEKEYEKAKKTYMQACKRSPSCLTWLGLGIACYRLEELTEAEDALSEANALNNYNAEVWAYLALVCLKVGRQLEAEQAYKYMILKCLKDEALLAEIHT
LQETVGFGNPSF

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

TTC18 (CFAP70) (NM_145170) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.