

Human Annexin A8-like Recombinant Protein 1/ANXA8L2 Recombinant Protein

Catalog Number: PROTQ5VT79

Overview

Product Name	Human Annexin A8-like Recombinant Protein 1/ANXA8L2 Recombinant Protein
Description	Human Annexin A8-like Recombinant Protein 1/ANXA8L2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-327aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q5VT79, 1-327aa; MGSSHHHHHH SSSLVPRGSH MGSMAWWKA WIEQEGVTVK SSSHFNPDPA AETLYKAMKG IGTNEQAIID
VLTKRSNTQR QQIAKSFKAQ FGKDLTETLK SELSGKFERL IVALMYPPYR YEAKELHDAM KGLGTKEGVI IEILASRTKN QLREIMKAYE
EDYGSSLEED IQADTSGYLE RILVCLLQGS RDDVSSFVDP ALALQDAQDL YAAGENIRGT DEMKFITILC TRSATHLLRV FEEYEKIANK SIEDSIKSET
HGSLEEAMLT VVKCTQNLHS YFAERLYYAM KGAGTRDGTI IRNIVSRSEI DLNLIKCHFV KMYGKTLSSM IMEDTSGDYK NALLSLVGSD P

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Annexin A8-like Recombinant Protein 1/ANXA8L2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.