

MDH1 Malate Dehydrogenase Chicken Recombinant Protein

Catalog Number: PROTQ5ZME2

Overview

Product Name	MDH1 Malate Dehydrogenase Chicken Recombinant Protein
Description	The DNA encoding Malate (Malic) Dehydrogenase is cloned from cDNA library of chicken heart. The MDH1 is purified by proprietary chromatographic techniques.
Size	1000IU, 5000IU, 50000IU
Tag	
Form	Sterile lyophilized powder.
Source	Escherichia coli
Formulation	Each mg of protein contains 0.59mg NaPO4.

Concentration

Storage

Lyophilized Malate dehydrogenase although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution MDH1 should be stored at 4°C between 2-7 days and for future use below -18°C. Please prevent freeze-thaw cycles.

Purity

Greater than 95.0% as determined by (a) Analysis by RP-HPLC and (b) Analysis by SDS-PAGE.

!description_after_attributes! ?/description_after_attributes?

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



MDH1 Malate Dehydrogenase Chicken Recombinant Protein