

Human B7-H3 (4Ig)/B7-H3b Recombinant Protein

Catalog Number: PROTQ5ZPR3

Overview

Product Name	Human B7-H3 (4Ig)/B7-H3b Recombinant Protein
Description	Human B7-H3 (4Ig)/B7-H3b Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 29-466aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q5ZPR3, 29-466aa; LEVQVPEDPV VALVGTDATL CCSFSPEPGF SLAQLNLIWQ LTDTKQLVHS FAEGDQGS A YANRTALFPD
LLAQGNASLR LQRVRADEG SFTCFVSIRD FGSAAVSLQV AAPYSKPSMT LEPNKDLRPG DTVTITCSSY QGYPEAEVFW QDGQGVPLTG
NVTTSQMANE QGLFDVHSIL RVVLGANGTY SCLVRNPVLQ QDAHSSVTIT PQRSPTGAVE VQVPEDPVVA LVGTDATLRC SFSPEPGFSL
AQLNLIWQLT DTKQLVHSFT EGRDQGSAYA NRTALFPDLL AQGNASRLQ RVRVADEGSF TCFVSIRDFG SAAVSLQVAA PYSKPSMTLE
PNKDLRPGDT VTITCSSYRG YPEAEVFWQD GQGVPLTGNV TTSQMANEQG LFDVHSVLRV VLGANGTYSC LVRNPVLQQD AHGSVTITGQ
PMTFPPEALE HHHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human B7-H3 (4Ig)/B7-H3b Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.