

Mouse Mer Recombinant Protein

Catalog Number: PROTQ60805

Overview

Product Name	Mouse Mer Recombinant Protein
Description	Mouse Mer Recombinant Protein expressed in Baculovirus with hlgG-His-tag. Sequence domain: 19-497aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, Q60805, 19-497aa; GGTAEKWEET ELDQLFSGPL PGRLPVNHRP FSAPHSSRDQ LPPPQTGRSH PAHTAAPQVT STASKLLPPV AFNHTIGHIV LSEHKNVKFN CSINIPNTYQ ETAGISWWKD GKELLAGHHS ITQFYDDEEG VSIIALFSIA SVQRSDNGSY FCKMKVNNRE IVSDPIYVEV QGLPYFIKQP ESVNVTRNTA FNLTCQAVGP PEPVNIFWVQ NSSRVNEKPE RSPSVLTVPG LTETAVFSCE AHNDKGLTVS KGVHINIKVI PSPPTEVHIL NSTAHSILVS WVPGFDGYSP LQNCISIQVKE ADRLSNGSVM VFNTSASP HL YEIQQQLQALA NYSIAVSCRN EIGWSAVSPW ILASTTEGAP SVAPLNITVF LNESNNILDI RWTKPPIKRQ DGELVGYRIS HVWESAGTYK ELSEEV SQNG SWAQIPVQIH NATCTVRIAA ITKGGIGPFS EPVNIHPEH SKVDYAPSST PPGNTDSM

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Mer Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.