

Mouse Serpin F2/alpha 2-Antiplasmin Recombinant Protein

Catalog Number: PROTQ61247

Overview

Product Name	Mouse Serpin F2/alpha 2-Antiplasmin Recombinant Protein
Description	Mouse Serpin F2/alpha 2-Antiplasmin Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 28-491aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

1mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, Q61247, 28-491aa; VDLPGQQPVS EQAQQKLPLP ALFKLDNQDF GDHATLKRSP GHCKSVPTAE ETRRLAQAMM AFTTDLFSLV
AQTSTSSNLV LSPLSVALAL SHLALGAQNG TLHSLHRVLH MNTGSLPHL LSHFYQNLGP GTIRLAARIY LQKGFPIKDD FLEQSERLFG
AKPVKLTGKQ EEDLANINQW VKEATEGKIE DFLSELPDST VLLLLNAIHF HGFWRTKFDP SLTQKDFHHL DERFTVSVDM MHA VSYPLRW
FLLEQPEIQV AHFPFKNNMS FVVVMPTYFE WNVSEVLANL TWD TLYHPSL QERPTKVWLP KLHLQQQLDL VATLSQLGLQ ELFQGPDLRG
ISEQNLVVSS VQHQSTMELS EAGVEAAAAT SVAMNRM SLS SFTVNRPLFL FIMEDTIGVP LFGVSVRNP N PSALPQLQE Q RDSPDNRLIG
QNDKADFHGG KTFGPD LKLA PRMEEDYPQF SSPKHHHHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Serpin F2/alpha 2-Antiplasmin Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.