

Mouse 5'-Nucleotidase/CD73 Recombinant Protein

Catalog Number: PROTQ61503

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Mouse 5'-Nucleotidase/CD73 Recombinant Protein
Description	Mouse 5'-Nucleotidase/CD73 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 29-551aa. Application(s): Enzyme Activity, SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, Q61503, 29-551aa; WELTILHTND VHSRLEQTS DSTKCLNASL CVGGVARLFT KVQQIRKEEP NVFLDAGDQ YQGTIWFTVY KGLEVAHFMN ILGYDAMALG NHEFDNGVEG LIDPLLRNVK FPILSANIKA RGPLAHQISG LFLPSKVLSV GGEVVGIVGY TSKETPFLSN PGTNLVFEDE ISALQPEVDK LKTLNVNKII ALGHSGFEMD KLIAQKVRGV DIVVGGHSNT FLYTGNPPSK EVPAGKYPI VTADDGRQVP VVQAYAFGKY LGYLKVEFDD KGNVITSYGN PILLNSSIPE DATIKADINQ WRIKLDNYST QELGRTIVYL DGSTQTCRFR ECNMGNLICD

AMINNNLRHP DEMFWNHVSM CIVNGGGIRS PIDEKNNGTI TWENLAAVLP FGGTFDLVQL KGSTLKKAFE HSVHRYGQST GEFLQVGGIH
VVYDINRKPW NRVVQLEVLC TKCRVPIYEP LEMDKVYKVT LPSYLANGGD GFQMIKDELL KHDSGDQDIS VVSEYISKMK VVYPAVEGRI
KFSLEHHHHH H

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Mouse 5'-Nucleotidase/CD73 Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.