

Mouse Sonic Hedgehog/Shh Recombinant Protein

Catalog Number: PROTQ62226

Overview

Product Name	Mouse Sonic Hedgehog/Shh Recombinant Protein
Description	Mouse Sonic Hedgehog/Shh Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 25-198aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, Q62226, 25-198aa; MCGPGRGFGK RRHPKLTPL AYKQFIPNVA EKT LGASGRY EGKITRNSER FKELTPNYNP DIIFKDEENT
GADRLMTQRC KDKLNALAI S VMNQWPGVKL RVTEGWDEDG HHSEESLHYE GRAVDITSD RDRSKYGMLA RLAVEAGFDW VYYESKAHIH
CSVKAEN SVA AKSGGLEHHH HHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Sonic Hedgehog/Shh Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.