

Rat GFR alpha-1 Recombinant Protein

Catalog Number: PROTQ62997

Overview

Product Name	Rat GFR alpha-1 Recombinant Protein
Description	Rat GFR alpha-1 Recombinant Protein expressed in Baculovirus with hlgG-His-tag. Sequence domain: 25-430aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Rat, Q62997, 25-430aa; DRLDCVKASD QCLKEQSCST KYRTLRCQVA GKETNFSLTS GLEAKDECERS AMEALKQKSL YNCRCKRGMK
KEKNCLRIYW SMYQSLQGND LLEDSPYEPV NSRLSDIFRA VPFISDVFQQ VEHISKGNNC LDKKACNLD DTCKKYRSAY ITPCTTSMNS
EVCNRRKCHK ALRQFFDKVP AKHSYGMLFC SCRDIACER RRQTIVPVCS YEERERPNCL SLQDCKTNY ICRSLADFF TNCQPESRSV
SNCLKENYAD CLLAYSGLIG TVMTPNYVDS SLSVAPWCD CSNSGNDLED CLKFLNFFKD NTCLKNAIQA FGNGSDVTMW QPAPPVQTTT
ATTTTAFRVK NKPLPGAGSE NEIPHTVLPV CANLQAQKLL SNVSGSTHLC LSDSDFGKDG LAGASSLEPK SCDKTHTCPP CPAPELLGGP
SVFLFPPKPK DTLMISRTPE VTCVVVDVSH EDPEVKFNWY VDGVEVHNAK TKPREEQYNS TYRVVSVLTV LHQDWLNGKE YKCKVSNKAL
PAPIEKTISK AKGQPREPQV YTLPPSRDEL TKNQVSLTCL VKGFYPSDIA VEWESNGQPE NNYKTPPVV DSDGSFFLYS KLTVDKSRWQ
QGNVFCSSVM HEALTHHNYTQ KSLSLSPGKH HHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Rat GFR alpha-1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.