

Human recombinant PLCXD3 protein

Catalog Number: PROTQ63HM9-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant PLCXD3 protein
Description	Human recombinant PLCXD3 protein expressed in E.coli with His-Tag. Sequence domain: 1-321aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E.coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q63HM9, 1-321aa; MGSSHHHHHH SGLVPRGSH MGSMASQK NELKLADWMA TLPESMHSIP LTNLAIPGSH DSFSFYIDEA
SPVGPEQPET VQNFVSVFGT VAKKLMRKWL ATQTMNFTGQ LGAGIRYFDL RISTKPRDPD NELYFAHGLF SAKVNEGLEE INAFITDHHK
EVVFLDFNHF YGMQKYHHEK LVQMLKDIYG NKMCPAIFAQ EVSLKYLWEK DYQVLVYHS PVALEVPFLW PGQMMPAPWA NTTDPEKLIQ
FLQASITERR KKGFFISQV VLTPKASTVV KGVASGLRET ITERALPAMM QWVRTQKPGE SGINIVTADF VELGDFISTV IKLNYVFDEG EANT

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant PLCXD3 protein](#)

For Research Use Only. Not for use in diagnostic procedures.