

Rat IL-33 Recombinant Protein

Catalog Number: PROTQ66H70

Overview

| | |
|--------------|--|
| Product Name | Rat IL-33 Recombinant Protein |
| Description | Rat IL-33 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 109-264aa. Application(s): SDS-PAGE. |
| Size | Starting from 20ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol |

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Rat, Q66H70, 109-264aa; MGSSHHHHHH SGLVPRGSH MGSSIQGTSL LTESCALSTY NDQSVSFVLE NGCYVINVED CGKNQEKDKV
LLRYESSFP AQSGDGVLDGK KLMVNMSPDK DTDIWLNAND KDYSVELQKG DVSPPDQAFF VLHKKSSDFV SFECKNLPGT YIGVKDNQLA
LVEENDESCN NIMFKLSKM

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Rat IL-33 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.