

B7 H6 (NCR3LG1) (NM_001202439) Human Recombinant Protein

Catalog Number: PROTQ68D85

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	B7 H6 (NCR3LG1) (NM_001202439) Human Recombinant Protein
Description	Purified recombinant protein of Homo sapiens natural killer cell cytotoxicity receptor 3 ligand 1 (NCR3LG1), with C-terminal His tag, secretory expressed in HEK293 cells, 50ug
Size	50 µg
Tag	C-His
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	PBS, pH7.4, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

A DNA sequence encoding the region Asp25-Ser262 of NCR3LG1

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

B7 H6 (NCR3LG1) (NM_001202439) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.